


CE 0459



PREVOR

ANTICIPATE AND SAVE
Toxicology Laboratory & Chemical Risk Management



The chemical burn




■ The aggressive chemicals

Two categories of chemicals can provoke a burn by a direct contact: corrosives and irritants. The seriousness of the chemical burn depends on the number of cells destroyed and will be more significant in the case of corrosives. In some cases, the aggressive product can also have toxic or harmful effects, the consequences of which should not be neglected.

CORROSIVE



Concentrated Acids and bases



SIGNIFICANT REACTIONS



IRREVERSIBLE EFFECTS

IRRITANT



Solvents, oils...



WEAK REACTIONS



REVERSIBLE EFFECTS




Risk of generalised toxicity

■ The mechanism of a chemical burn

THE EXCHANGE


SKIN and EYE

CHEMICALS



A chemical burn is caused by the capability of corrosives and irritants to create an exchange (ion, proton, electron...) with the tissue of the skin or the eye. The degree of the burn will depend on the number of molecules destroyed and on the type of modification (reversible or irreversible)

■ The action of soda NaOH:



Healthy cells as seen under a microscope

Addition of soda 0.5N (2%)

Cells are completely destroyed

* Picture source: ACTO, Prof. Norbert Schrage, Aachen.

The principles of emergency washing



To stop the aggressive reactions

- There are 6 types of aggressive chemical reactions



• Acid •



• Basic •



• Oxidising •



• Reducing •



• Chelating •

ACETONE
ETHANOL

• Solvent •


To be effective regardless of the type of accident and to avoid the risk of error for the victim, it is necessary to be able to stop these 6 aggressive reactions.



A POLYVALENT PRODUCT IS ESSENTIAL



To stop the evolution of the burn



A chemical burn is initiated by the contact between the aggressive product and the skin or eye. Following this contact, part of the aggressive chemical will penetrate into the tissue and cause destruction of the cells.

To effectively decontaminate the splash, it will be necessary, not only to decontaminate the product on the surface, but also to control its penetration inside the tissue.



A PRODUCT ABLE TO STOP THE AGGRESSIVE CHEMICAL IS ESSENTIAL

■ Factors influencing the penetration

- The type of product and its concentration
- The temperature
- The length of time of exposure




H₂O




First-aid treatment: From water...




- **The principles of washing with water**
 - washing of the surface to remove quickly the aggressive product
 - dilution of the chemical to reduce its aggressiveness
 - a universal product avoids the risks of error at the time of the accident
- **What are its limitations?**
 - concentrated products which penetrate very quickly
 - the intervention time of 10 seconds which is not always achievable
 - washing comfort : risk of hypothermia under a water shower
difficulty to open the eye



Healthy cells as seen under a microscope



The beginning of washing: water penetrates inside the cells and makes them expand.



The end of washing : cells are destroyed because water causes them to burst.



In vitro tests

The effectiveness of Diphoterine® has been tested on more than 600 chemicals belonging to the various classes of aggressive products (acids, bases, oxidising agents, etc...) within the Prevor laboratories. This list is available on our web site.



In vivo tests

Application of soda on rabbits' corneas:



The ideal product in case of

Retain the advantages of water:

- Fast washing of the surface
- Single protocol

Washing with water: 100% of the endothelial cells are destroyed

Washing with the Prevor solution : 95% of the endothelial cells are preserved

P. Josset, B. Belosse, J. Blomet, H. Saraux, Intérêt d'une solution isotonique dans le traitement précoce des brûlures chimiques basiques cornéo conjonctivales. Bull Soc Opt France 1986, 6-7, LXXXVI.



...to Diphoterine®


■ The principles of washing with Diphoterine®

- It is a liquid which enables it to obtain the same effect as water on the surface of the skin or eye
- Diphoterine® is a chelating amphoteric agent, which enables it to stop the aggressiveness of the chemicals in a polyvalent way. *
- The Diphoterine® mechanism can be illustrated as follows:




BASE ACID

Diphoterine® will attract the chemical in contact with the tissue



BASE ACID

The acid site of Diphoterine® will fix base to make them harmless



BASE ACID


The basic site of Diphoterine® will fix acid to make them harmless

- Diphoterine®, unlike water, will **stop the penetration** of the chemical.


After a chemical splash must:

Bring improvements to compensate for the limitations of water:


- guarantee a total efficacy whatever the product
- Increase the intervention time
- Improve washing comfort to increase the effectiveness



Healthy cells as seen under a microscope.



Beginning of the washing with Diphoterine®: the cells contract slightly




End of the washing: the cells are preserved.

Results

Series of cases in industry:

The introduction of Diphoterine® enabled a reduction in the average duration of days lost from work to a few hours. The standard deviation experienced in this case reduced from 8.12 to 0.4, which illustrated that Diphoterine® provided a more reliable first-aid treatment following a chemical splash.



The effectiveness of Diphoterine® has been proven in industry through many reported and published cases of accidents (complete bibliography available on Internet) At Martinswerk (Germany), the effectiveness of Diphoterine® was compared with the effectiveness of water on 45 splashes with soda. The results can be obtained in the following table:



How to use

A sterile washing solution

Dispensers for the skin, in the form of micronised sprays (optimising the surface of contact)



> To decontaminate an entire body

Typical equipment for locations where there are large quantities of chemicals :

- Production areas
- Warehousing areas
- Off-loading and decanting areas

> To decontaminate a limb or a torso

An ideal complement for a Wall Mounted Eyewash, where the quantities of chemicals are limited:
Laboratories

> To decontaminate a face or a hand

Ideal personal equipment for small splashes:
Workshops
Maintenance personnel



Washing must begin within the first minute of the splash

Protocol for the first aid in the event of a splash



1/ Go away from the danger



2/ Get undressed



The Diphoterine[®] ?

To be used immediately after the accident

Dispensers for the eyes, equipped with an ergonomic eyecup, which assists the opening of the eye for a more effective washing

LPM
Portable
mural
eyewash



Portable
Eyewash



SIEW



> To decontaminate an eye

In areas where the equipment must be portable:

- First-aid boxes
- Medical centre
- Emergency vehicles



> Personal equipment to decontaminate an eye:

Ideal for the maintenance personnel



> To decontaminate an eye and the skin

In locations where fixed equipment is well adapted:

- Laboratories
- Production areas
- Warehousing areas
- Off-loading and decanting areas



Washing must begin within the first 10 seconds of the splash

Washing must begin within the first minute of the splash

The first aid intervention with Diphoterine[®] in the event of a chemical splash



3/ Wash as quickly as possible, respecting the instructions for use of Diphoterine[®]



4/ Alert



5/ Seek medical advice






PREVOR

Science to improve the safety

Our Missions

> 1/To understand the chemical products

Research on the toxicity of the products
Studies on the mechanism of the burn



> 2/To find solutions

- => Diphotérine®
- => Hexafluorine® (HF)
- => Trivorex®

> 3/To advise

- Our customers for a more effective use (advice on site, training of the users and consultants, edition of informative works regarding our knowledge)
- Institutions involved with the chemical risk

Reference organisations They trust us

3M FRANCE - ADP - AIR FRANCE - AKZOL NOBEL - ALCATEL - ARJO WIGGINGS - ASTRA - ATOCHEM - BASF PEINTURES ET ENCREs - BAYER POLYSAR FRANCE - BEGHIN SAY - BOEHRINGER INGELHEIM - BP CHEMICALS - CALAIRE CHIMIE - CALCIA - CEA - CITERNORD - COCA COLA PROD - COGEMA - COMPAGNIE GÉNÉRALE DES EAUX - CRISTAL UNION - CRISTALLERIE D'ARQUES - DANONE - DE CARBON - DE DIETRICH - DUPONT DE NEMOURS - EDF - ELF - ESSILOR INTERNATIONAL - GEC ALSTHOM - GENERAL ELECTRIC - GENERAL DES EAUX GLAENZER SPICER - GLAXO - GOODYEAR - GSK - GUERBET - GUERLAIN

HEINEKEN - HENKEL FRANCE - HBL - HOËCHST L'AIGLE - ICI - JAGUAR CARS LABORATOIRES ECLAIR - LIPHA - L'ORÉAL - LUBRIZOL FRANCE - LYONNAISE DES EAUX - MICHELIN - MUSÉE D'ORSAY-LOUVRE - ORIL-SERVIER - PASTEUR VACCINS PEUGEOT - PFIZER - PHILIPS CI - PPG INDUSTRIES FRANCE - RECKITT ET COLLMAN RENAULT - RHÔNE POULENC - ROHM AND HAAS FRANCE - ROQUETTE FRÈRES SAFT - SAINT-GOBAIN - SCHLUMBERGER - SNCF - SOLLAC - SOLWAY ET CIE SOPAD NESTLÉ - TRW - UGINE - VALLOUREC...



www.prevor.com

Moulin de Verville
95760 Valmondois
Tel. 00 33 (0)1 30 34 76 76
Fax. 00 33 (0)1 30 34 76 70